# Reactivos GPL

Barcelona, España

-LABTROL H NORMAL-  
-LABTROL H  
PATHOLOGICAL-

CONTROL "H" NORMAL  
/NORMAL  
CONTROL "H" PATHOLOGICAL  
/PATOLÓGICO

Store at : +2+8°C.  
Conservar entre: +2+8°C.

Presentation / Presentación:  
Cod. QC003 Control Human Normal 4 x 5  
Cod. QC004 Control Human Pathological 4 x 5

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**Serum control multicomponent for clinical chemistry assays.**

**Only for in vitro use in clinical laboratory (IVD)**

### CHARACTERISTICS OF PRODUCT

Human lyophilized serum.  
With most constituent concentrations and activities in the normal or pathologic range.

IS INTENDED FOR CONTROL OF ACCURACY FOR USE WITH  
MANUAL AND AUTOMATED ANALYTICAL PROCEDURES.

### REAGENTS COMPOSITION

Human serum. Enzymes (human, bovine and porcine). Nonprotein constituents and bacteriostatic agents. The concentration / activities of the components are lot-specific. The exact values and ranges valid for reagents are given informational purpose in the value sheet.

### PRECAUTIONS

The individual units comprising lot or lots of components from human origin, have been tested for HBsAg and antibodies to HIV and HCV found negative / non-reactive by an FDA approved test. However, handle cautiously as potentially infectious. Avoid contact with skin and eyes. In case of contact with eyes, rinse immediately with plenty of water and seek medical device. Wear suitable gloves and eye/face protection. In case of an accident, seek medical device immediately and show this container or label.

### REAGENT PREPARATION AND STABILITY

- Using Class A volumetric pipette, carefully add the exact amount of distilled or deionized water as specified on the vial label to the control material.
- Swirl the vial gently and allow to stand at room temperature for 5 to 10 minutes.
- Swirl content until homogeneous. DO NOT VORTEX. Remove the sample required for testing.

Improper handling and/or storage of the control can affect results. Inaccurate reconstitution of the control and errors in assay technique can cause erroneous results.

The Control serum is stable until the expiration date on the label when stored tightly closed at 2-8°C and contaminations are prevented during their use.

Do not use reagents over the expiration date or if there is visible evidence of microbial growth. Store tightly capped when not use.

- After reconstitution is stable for:  
At 15°C to 25°C: 12 hours (\*8 hours Bilirubin Total; 4 hours Bilirubin Direct, ACP and ALP).  
At 2-8°C: 5 days (\*24 hours Bilirubin Total, Bilirubin Direct, ACP and ALT).  
At -15°C to -25°C: 1 month (\*2 weeks Bilirubin Total, Bilirubin Direct, ACP and ALP).
- Alkaline phosphatase may increase in activity with time.

### BIBLIOGRAPHY

1. Council Directive (2000/54EC). Official Journal of the European Communities No. L262 from Oct. 17th, 2000.
2. International Federation of Clinical Chemistry (IFCC) Education Division, Expert Panel of Quantities and Units: A Protocol for the Conversion of Clinical Laboratory data, Journal of Automatic Chemistry Vol.11, No.5 (Sept-Oct 1989), pp.223-226

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CHEMELEX, S.A.